



Innovative, environmentally aware thinking, safety, and positive patient outcomes are central to all we do. Pintler Medical leads the industry in creating groundbreaking active warming and pressure redistribution technology.

We know that by improving the comfort and circulation of patients during the surgical procedure, we can dramatically improve the clinical outcomes as well as the patient's personal satisfaction.

Patient satisfaction delivers positive surveys which maximize your Active Warming Reimbursement.\*

"Finally, a warming system that actually makes our patients warm!! Whatever you're doing, don't stop!"

- Chris, OR Director, VA Hospital, Providence, RI

"Saves me steps and time. And is much better for the patient...I'm all for it!" - Janie, CRNA, Lawrence, KN

"When we performed a skin survey 24 hours **after** a 17 hour procedure, absolutely no skin issues were detectable. We were all pleasantly surprised!"

- RN, St Luke's, Cedar Rapids, IA

"All of our patients are warmer now than they have ever been."
- PACU, Palos Community Hospital, Palos Heights, IL

"We love these pads!"

- OR Staff, Sacred Heart, Spokane, WA

It was so much easier to turn our rooms over, we didn't have to prep or do any extra work."

- Durham VA, Durham, NC

\*CMMS.Gov: PQRS Measure 193

SCIP-10 v3.2 Resistive Warming is Active Warming

## PINTLER PEAK™ OR Warming Pads

Pintler Peak OR Warming Pads maintain patient normothermia and pressure redistribution throughout surgery. The active warming technology delivers positive warming immediately upon contact and maintains warming throughout the procedure and safely heats only the patient. Pressure redistribution encourages circulation and provides reduced risk of pressure ulcers.

- Resistive warming does not disrupt the downward flow of air in surgery
- Pintler operating controls easily mount to the anesthesia IV pole
- · Sized to fit any OR table model or measurement
- Waste-reducing, anti-microbial and antibacterial patient surface
- Minor punctures are repairable with Pintler Repair Kit
- Ease of use allows for utilization during all procedures, even those less than 1 hour
- Increased opportunity for Active Warming reimbursement
- Safe for use with Pediatrics

## PINTLER VISTA™ Heated Stretcher Pad

Pintler Vista Heated Stretcher Pads assist in providing a normothermic PreOp/PACU environment. An innovative pendent design gives the patient the ability control their temperature setting, reducing the use of blankets and eliminating disposable warming products. Patients personal comfort and consistent warmth translate into increased patient satisfaction and positive survey results.

- Patient comfort when a patient remembers it most:
   right before and right after a procedure
- Patient-controlled warming frees nursing staff to focus on patient care and safety
- Patient satisfaction increases positive surveys
- Fits every patient stretcher model and measurement
- Works well for clinical and diagnostic procedures



## PINTLER BASE+ Unheated Surgical Pads

Patient satisfaction increases positive surveys, which result in improved hospital reimbursement. From the moment a patient is transferred to a Pintler Base+ Pad, they can be assured of pressure relief. The construction of our pads allows for maximum weight disbursement.

- Our superior quality foam and OrcaSeal™ coating provide maximum pressure redistribution
- Minor punctures repairable with Pintler Repair Kit
- Pad foam and OrcaSeal™ coating are anti-microbial and anti-bacterial
- All pads are designed to meet all OR table configurations
- Bariatric and specialty pads available
- Unique pad design removes the need to use Velcro



## PINTLER MEDICAL, LLC

Seattle, WA/Philipsburg, MT

www.pintlermedical.com • info@pintlermedical.com

877.620.9868 • 206-268-0441

Veteran Owned Small Business
All Products Made in the USA

D . I	• 1		- 1	1	
Distr	ın	$\Box$	$\sim$	$n \cdot /$	•
レスロ	$\mathbf{I}$	$\alpha \iota c$	u	$\cup$ $\vee$	